




CLINICA

- ID:
- Name:


CASO

- ID:
- OPERADOR: ESCANER LPP
- FECHA :2018 16/05/2018


MODELO ORIGINAL




MOVIMIENTO 1




MOVIMIENTO 2



MOVIMIENTO 3




MOVIMIENTO 4





ANALISIS MOVIMIENTOS

Overview



ANALISIS DE MOVIMIENTOS

Setup: SUP

Tooth Number	Rotation	Angulation	Inclination	Left/ Right	Extrusion/ Intrusion	Forward/ Backward	IPR Distal (mm)	IPR Mesial (mm)	Long Axis	Circumference
12			5.2 deg	0.3 mm	0.0 mm	0.4 mm			86.6 grade	21.2 mm
11			1.3 deg	0.1 mm	0.0 mm	0.5 mm			86.6 grade	25.4 mm
21			1.6 deg	0.0 mm	0.0 mm	0.5 mm			88.6 grade	26.2 mm
22			1.0 deg	0.3 mm	0.0 mm	0.4 mm			84.8 grade	20.9 mm

Setup: SUP Subsetup1

Tooth Number	Rotation	Angulation	Inclination	Left/ Right	Extrusion/ Intrusion	Forward/ Backward	IPR Distal (mm)	IPR Mesial (mm)	Long Axis	Circumference
--------------	----------	------------	-------------	-------------	----------------------	-------------------	-----------------	-----------------	-----------	---------------